

NOTICE: This document contains references to Varian.  
Please note that Varian, Inc. is now part of Agilent  
Technologies. For more information, go to  
**[www.agilent.com/chem](http://www.agilent.com/chem)**.



## APPLICATION NOTE

### β-Blockers

Application 2326 LC

Technique : HPLC  
Column : Pursuit™ 5μ C8, 150 x 4.6mm  
Part no. : Conventional A3030150X046  
: ChromSep Cartridge System A3030150C046  
Temperature : Ambient  
Eluent : A: MeCN  
: B: 0.1% NH<sub>4</sub>OH in water, pH 10.3  
: A:B = 65:35  
Flow rate : 1.0 mL/min  
Detection : UV @ 215 nm  
:

Courtesy : Varian Inc.

### Peak Identification

1. Hydrochlorothiazide
2. Practolol
3. Pindolol
4. Bisoprolol
5. Alprenolol

